
**ENG** Portable Oral Irrigator AQUAJET LD-A3 with accessories  
Instruction Manual

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ENG			
BEFORE USING THE DEVICE, READ THIS OPERATION MANUAL CAREFULLY AND SAVE IT FOR FUTURE USE. Manual is described with the following terms:			
<b>Warning</b> Requirements, non-compliance with which may lead to injuries or property damage.			
<b>Caution</b> Forbidden actions and device handling methods.			
<b>Note!</b> Mandatory actions and instructions, which should be strictly followed to ensure safe operation.			
MAIN PARTS			
Figure №	Name	Figure №	Name
1.1	Nozzle (2 pcs, standard LD-SA01 nozzle)	2.5	Power cord plug
1.2	Nozzle disconnect button	3.1	Short tube
1.3	Mode switch	3.2	Short tube filter
1.4	Indicator	3.3	Container
1.5	Tube jack	4.1	Connector jack for long tube
1.6	Power supply jack	4.2	Connector
2.1	Power supply	4.3	Long tube
2.2	Power supply USB-port	4.4	Long tube holder
2.3	USB-plug for power cord	4.5	Long tube filter
2.4	Power cord		

The main body and nozzle are mainly made from ABS, the container is made from PP. The nozzle has passed the biocompatibility test. People who are sensitive to such materials should take care when using the device.

**INTENDED USE**

AQUAJET Portable oral irrigator is intended for treatment and prevention of periodontal diseases, gingivitis, periodontitis; for medical irrigation and for the massage of the oral cavity, and care for dental appliances (including crowns, implants, prosthetic dentures, occlusion correction devices) and interdental spaces.

**AQUAJET usage for oral hygiene cannot substitute the toothbrush completely.**

**Potential consumers:** all people older than 2 years.  
**Functional features:** using a powerful pulsating water-jet, the device removes food debris from interdental spaces and other difficult areas in the oral cavity, massages gums, improves the blood flow, and prevents gum bleeding. It provides maximum care for prosthetic dentures and crowns, prevents the development of pathogenic bacteria and inflammatory processes in places of prosthetic denture and crown interface with the oral mucosa, which is not achieved by other common means. It prevents the formation of dental calculus, gingivitis, and many other periodontal diseases. Compact size and power supply from a self-contained battery allow the usage of the irrigator while travelling. The battery is charged from the power supply (PS) or another source equipped with a USB-port and corresponding to the characteristics defined in the section "Technical characteristics".

**Use of AQUAJET Irrigator by Children**  
**Device usage by a child is allowed only under the surveillance of an adult who has thoroughly read this Operation Manual.** Irrigator usage by children decreases the risk of dental diseases, facilitates oral care, and shortens oral hygiene time significantly in case of brace wearing. Children also often lack patience and diligence for thorough tooth cleaning, while the irrigator usage provides the best oral hygiene possible. Children may start using oral irrigator when they can use the toothbrush independently. Though correct usage of the oral irrigator cannot traumatize the oral cavity, pay a visit to the pediatric dentist before using the irrigator.

**Device usage by a child is allowed only under the surveillance of an adult who has thoroughly read this Operation Manual.** Infants and children without supervision should stay away from the device, as they might swallow the small parts of the device. Supervise children and infants carefully while the device is working.

**SAFETY MEASURES**

For safe device usage, carefully read this Operation Manual. Non-compliance with the section requirements may lead to hazards of various grades for the user and people around the user.

When using electric appliances, especially in the presence of children, always follow the general precautions mentioned below.  
Do not handle the plug with wet hands. Before charging the device, check that the power supply cord is not damaged. This device is equipped with a built-in battery. Do not throw the device into the fire, do not heat it, do not use or leave it in rooms with high temperature.

Do not place the device into the water. Do not use it when bathing. Do not immerse the power supply into the water and do not wash it with water. The power supply cord should not touch hot or heating surfaces. If the device is not working, do not try to repair it. Contact authorized repair services. Do not leave the switched-on device unattended. Use only the appliances which are intended for usage with this device and are described in this Operation Manual. Do not use accessories which are not recommended by the manufacturer. Do not insert foreign objects into the device holes. Use the device only for the oral care. Water head applied to eyes, nose, ears, or throat may lead to injuries. Do not use water with temperature exceeding 40°C for the device. Do not service and maintain the device while it is being used by the patient. Keep it away from the pets and pests to avoid the contamination by them. This device don't used with HF (High Frequency) surgical equipment at the same time. To avoid problem with the device, do not use the device near strong electromagnetic fields, radiated interference signal or electrical fast transient/burst signal. For example magnets, radio transmitters, microwave ovens. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm to any part of the dental oral irrigator, including cables specified by the manufacturer.

When disconnecting the power supply from the outlet, always hold the plug, not the cord. After using the device, empty the container with leftover water and turn on the device for 2-3 sec. for complete drying.  
**Risks of Usage**  
This device is not intended for usage by people with limited physical, sensory, or mental abilities, or persons without required knowledge and experience. The device usage by such people is allowed only under the surveillance of the person responsible for their safety, with him providing clear instructions on the device operation. Make sure that children do not play with the device.

**Contraindications:** patients with periodontal abscesses and oral ulcers, persons who cannot operate the device independently.  
**Expected side effects:** no expected side effects detected.

**BATTERY CHARGING**

The device uses a battery intended for multiple charging and prolonged use. Before the first device use, the battery should be charged for 16 hours. During charging, the main unit slightly heats, but it is not a defect. Charge the battery at room temperature within the 15-35°C range. When charging the battery at a lower temperature, the device operating time may decrease, and at temperatures below 0°C charging may become impossible. The manufacturer recommends using the PS (power supply) LD-N064 (included). Also, the included PS cord may be used for charge the device from any PS with a USB-port (type A). Recommended parameters are described in the "Technical characteristics" section.

To charge the battery, insert one end of the PS cord (2.3) into the port on the device back-panel (1.6), and the other end of the cord (2.3) — into the PS USB-port (2.2). Insert the PS plug into the wall outlet (Fig. 5). The charging process is indicated by the blue indicator on the device (1.4). When the battery is fully charged, the indicator turns to green.

For electric safety, the device cannot operate and charge the battery simultaneously.

The battery is charged fully approximately within 4-8 hours. At full charge, the device may operate approximately 20 minutes on maximum power. Low battery charge status is indicated by a yellow light.

**OPERATION PROCEDURE**

The use of alcohol-containing liquids, liquids with suspensions, oily fluids, and other fluids not intended for usage in the irrigator may lead to irrigator failure!

The manufacturer recommends using filtered tap water for irrigation. Internal parts of the new device may contain residues of distilled water, which was used for the device testing. This is not a defect and it is completely safe.

**Nozzle Connection/Disconnection**  
The set includes 2 individual nozzles with colored rims. Choose the individual nozzle of a certain color for each family member. Use only your nozzle. Insert the nozzle into the device until it engages, so that the nozzle's guide aligns with the groove on the irrigator handle (Fig. 6). To disconnect the nozzle, press the disconnection button and pull the nozzle, holding it in the notch area (Fig. 7). Recommendations for safe nozzle disconnection:

- If the nozzle disconnects from the irrigator poorly or slides in your hand because of water, use a dry cloth. Wrap the nozzle with a non-rough dry cloth and repeat the disconnection procedure.
- When disconnecting the nozzle from the irrigator, pull it WITHOUT rotation. Nozzle rotation may lead to its damage.

Additional nozzles may be purchased separately. Purchase only AQUAJET LD-SA01, AQUAJET LD-SA02 or AQUAJET LD-SA05 nozzles for this device. The irrigator may be used with an included container or another water vessel.

**Irrigator Use with the Included Container (Fig. 14)**  
Remove the connector (4.2). To do this, turn it clockwise, so that the indicator ▲ is in the "open" position (Fig. 8). Insert the short tube plug (3.1) into the jack in the lower part of the device (1.5) (Fig. 9). Fill the container with warm water to the MAX level. Connect the container to the device (Fig. 10) and turn it counterclockwise, so that the indicator ▲ is in the "closed" position (Fig. 11). The device is ready for use. Completely filled container lasts for 40 seconds of operation.

Do not use the device without water. If the device is tilted too severely when using the container, the tube filter (4.5) may escape the water, and water supply may stop.

**Double Aqua Flow system – Irrigator Use with Other Vessels (Fig. 15)**  
Double Aqua Flow system (DAF) allows usage of any available storage basins or standard container. If the container volume is not sufficient for you, the device may be used with other water vessels, e.g. a simple glass. To do this, insert the long tube plug into the jack (4.1) on the connector back side (4.2), if the connector is attached to the device (Fig. 11), or into the jack (1.5) in the lower part of the device, if the connector is disconnected from the device (Fig. 9). Unwind the long tube (4.3) and put the tube holder (4.4) onto the vessel with warm water, so that the tube filter (4.5) is just off the vessel's bottom (Fig. 12). The device is ready for use. Recommended continuous operation time should not exceed 5 min.

During the first use, pour warm water into the container, put the mode switch into position ●●●, so that the water fills all the device. The device will be ready for the first use when the water in the container is spent completely. Check that the temperature of the fluid in the container does not exceed 40°C. Do not use hot fluids, as the oral mucosa is very sensitive and may be damaged by a hot solution.

**Use**  
Bend over the sink. Bring the nozzle tip to the teeth or gums. Keep your mouth slightly opened, so that the water can pour into the sink. Direct the nozzle perpendicularly to farther teeth of the lower dental arch. Turn on the device, switching the mode switch (1.3) to position ●●● (low pressure) or to position ●●●● (normal pressure). Slowly and carefully, clean the interdental space for 4-5 sec., holding the jet near your gums. Then swipe the jet along the gum line near the tooth base to the next interdental space. Repeat the same movements as for the first interdental space (Fig. 16: dashed arrow indicates the place of water jet application, solid arrow — the nozzle movement direction).

You can stop the jet at any time, switching the mode switch into the position ○●● or ○●●●. Do not direct the jet directly into the gum pocket; try to hold the jet perpendicularly to the tooth line. Perform the same procedure for all interdental spaces from external and internal sides of the lower tooth arch. Then repeat the procedure for external and internal sides of the upper tooth arch. (Fig. 17) If artificial orthopedic structures are present in the oral cavity, their places of contact with gums, teeth, or other dental appliances should be rinsed with a jet more thoroughly. It is recommended to consult the dentist on the features of existing rinsing channels in these structures.

After 10 minutes of continuous operation, the device automatically turns off, and the green indicator starts flashing on and off. Do not use the device for at least 10 minutes to ensure cooling down of internal mechanisms.

Always empty the water from the container after usage. This prevents bacteria from multiplying. Switch off the device by switching the mode switch (1.3) into the position ○●●, remove the nozzle, the container, and the tube. Rinse everything under warm tap water and wipe them with a dry cloth. For convenient storage and transportation, put the accessories into the included container and connect it to the main unit (Fig. 13).

**BATTERY USE**

Description: the device uses a rechargeable battery (RB) intended for prolonged use. The RB operation period: 3 years in case of correct use. Recommendations: during operation, discharge the batteries completely to prevent battery memory effect, which shortens the operating time when fully charged. It is recommended to leave batteries charged at 70-80% for long-term storage. RB maintenance: not required.

RB replacement: if necessary, RB replacement is performed only in authorized service centers. To access the RB, remove the mode switch, disconnect the connector, unscrew 3 screws in the lower part of the main unit (under stops), and pull out the internal part. RB type and model: GP 130AAHV2B6.

**WARRANTY**

The warranty for the device is 12 months after the purchase date. Warranties are provided in the warranty card when selling the device to the customer.

**HANDLING, STORAGE, REPAIR, AND DISPOSAL**

- It is recommended to clean the device after each use.
- Before cleaning the device, disconnect the PS cord. Regularly clean the device with a moist cloth. The container, tubes, and nozzles may be washed in the dishwashing machine. Do not use rough cloths, brushes or abrasive agents for device cleaning.
- Use only soft washing agents. Never use solvents, petrol, or alcohol, as they may lead to operation failures or damage to device parts.
- Prepare the device from direct sunlight and impacts.
- Do not store and do not use the device in close proximity to heating devices or open fire.
- Avoid device contamination.
- Avoid device contact with aggressive solutions.
- If necessary, repair the device only in authorized organizations.
- This device is equipped with nickel metal hydride batteries. Please, ensure that battery disposal is performed in the official receiving center, if one is present in your region.

Do not modify this equipment without authorization of the manufacturer. All major maintenance of the device must be performed by an authorized service of authorized representative or manufacturer!

**TECHNICAL CHARACTERISTICS**

Model	LD-A3
Power consumption	5V, 500mA
Battery capacity	2x1250 mA/h
Battery life	approximately 3 years*
Water pressure:	
pressure regulator in position ●●● (low)	approximately 500 kPa (5 kgs./sq.cm)
pressure regulator in position ●●●● (normal)	approximately 900 kPa (9 kgs./sq.cm)
Pressure regulation	2 positions (low, normal)
Water jet pulsation frequency	1500 pulses/min.
Automatic switch-off time in continuous operation	10 min
Recommended max. continuous operation time	5 min
Time required for device to cool	at least 10 min
Time of complete vessel volume usage	approximately 40 sec.
Package contents	main unit, connector, container, 2 nozzles AQUAJET LD-SA01, long tube with holder, short tube, power supply LD-N064, power supply cord, operation manual, warranty card, individual carton package. 160 ml.
Container capacity, minimum	
Device operation conditions:	
ambient air temperature	from 10°C to 35°C
humidity	15% - 80% Rh
atmospheric pressure	from 86 to 106 kPa
Device storage and transport conditions:	
ambient air temperature	from minus 10°C to 40°C
humidity	15% - 95% Rh
atmospheric pressure	from 50 to 106 kPa
Noise level	no more 54** dB
Irrigator protection degree	IPX5
Power supply LD-N064 (included)	
Output voltage	5 V
Maximum load current	500 mA
Nominal voltage	100-240 v/50 Hz
Dimensions (main unit)	
with connector	215 x 50 x 54,5 mm
without container	297 x 50 x 54,5 mm
Weight (without package)	350 g
Year and month of manufacturing	indicated in the bottom of the unit body in a serial number after symbol "A" (first two digits - year, second two digits - month).

\* Battery life depends on the operation and storage conditions.  
\*\* Measured by Little Doctor International (S) Pte. Ltd.  
When the ambient temperature is 20°C, keep the device at room temperature for more than 2 hours before usage when the device is moved from low temperature storage environment or high temperature. The revision date of this User Manual is on the last page as XX-XXXX-YYMM-NN, where YY is the year, MM is the month, and NN is the revision number.  
Device manufacturing conforms to ISO 13485 international standard. The device is compliant with European Directive MDD 93/42/EEC, requirements of IEC 601-1-88, IEC 60601-1-2.

**Symbols:**  
Symbol for marking of electrical and electronics devices according to European Directive MDD 93/42/EEC  
IPX5 protection class  
Type B Applied part  
Electric shock protection class II  
Important: Read the manual  
When disposing of the device, consult your current regional regulations.  
Manufacturer  
Storage, transport and operation conditions  
Serial number

Indicator Lights:	Meaning
Light color	
Green light	The battery is fully charged. Run normally
Green flashing	Automatic switching off. It takes 10 minutes to cool the device.
Yellow light	Low battery charge
Blue light	During charging process


**TYPICAL FAILURES**

Problem	Cause	Solution
The device does not work	Battery is discharged	Charge the battery
	No water in the container	Fill the container with water
	Tube end is above the water level	Pour additional water so that the tube end is in water
	The tilt angle of the device is too big	Do not slope the device excessively
	The tube inserted into the socket not well.	Insert the tube into the socket well.
	Short/long tube's filter is clogged	Rinse the filter with tap water
	Long tube is twisted	Straighten the long tube

If, despite the aforementioned recommendations, you cannot achieve device operability, stop its operation and contact the organization performing maintenance (addresses and telephones of authorized organizations are given in the warranty card). Do not try to repair the irrigator by yourself.

**INFORMATION ON THE MANUFACTURER AND DISTRIBUTORS**  
Manufactured under control and for Little Doctor International (S) Pte. Ltd. (7500A BEACH ROAD #11-313 THE PLAZA SINGAPORE 199591. Postal address: Yishun Central P.O. Box 9293 Singapore 917699).  
Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No. 8, Tongxing Road Economic & Technical Development Area, 226010 Jiangsu, Jiangsu, People's Republic of China).  
Distributor in European Union: Little Doctor Europe Sp. z o. o. (S7G Zaplawa Street, 30-390 Kraków, Poland, phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53, e-mail: biuro@littledoctor.pl)  
For more information please visit www.aquajet.sg

**ADDITIONAL IRRIGATOR ACCESSORIES\***  
**AQUAJET LD-SA01 Nozzle**  
● Intended for use with AQUAJET LD-A3, LD-A8 irrigator  
● Made of plastic  
● For individual use  
● Number in the package — 2 pcs.  
**AQUAJET LD-SA02 Nozzle**  
● Special subgingival nozzle with a flexible nozzle for cleaning and hygiene of difficult areas  
● Intended for use with AQUAJET LD-A3, LD-A8 irrigator  
● Made of plastic with a flexible nozzle  
● For individual use  
● Number in the package — 2 pcs.  
**AQUAJET LD-SA05 Nozzle**  
● Intended for use with AQUAJET LD-A3, LD-A8 irrigator  
● Special orthodontic nozzle with a toothbrush for clearing dental brackets  
● Made of plastic  
● For individual use  
● Number in the package — 2 pcs.

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